

2-Propanol & Benzalkonium Chloride Antiseptic Solution

FUMISPRAY[®] - RAPID



FOR QUICK & EFFECTIVE DISINFECTION OF HIGH-TOUCH SURFACES, ELECTRONIC EQUIPMENT & SMALL OBJECTS



Composition:

Each 100 g Contains:
Isopropyl Alcohol IP 63.0 g
Benzalkonium Chloride 0.025 g

**Quick drying and effective
Kills > 99.99% germs in just 10 seconds**

Aldehyde free & safe for daily use - does not leave behind any residue or stain

Convenient pack to carry anywhere and spray

Ideal for keyboards, phones & other small objects & surfaces



**Packaging:
100 ml handy bottle with mist spray pump**



ISO 13485 Medical Devices
Certified QMS



PSK PHARMA PRIVATE LIMITED

121/22 Mahajenahalli, Shimoga Road, Harihar - 577601, Karnataka, India
Tel: +91 96060 69014 | e-mail: weprotect@pskpharma.com | www.pskpharma.com



WHO-GMP Certified

2-Propanol & Benzalkonium Chloride Antiseptic Solution **FUMISPRAY® - RAPID**



**FOR DISINFECTION OF SENSITIVE EQUIPMENTS &
SKIN PREPARATION**

**ALDEHYDE FREE!
DOES NOT CONTAIN ALDEHYDES**

Each 100 gms contains:

2-Propanol IP 63 Gms
Benzalkonium Chloride 0.025 Gms
(As Benzalkonium Chloride Solution IP)

Salient Features:

- Contains Isopropyl Alcohol 63% by weight
- Effective against bacteria, fungi, viruses
- Dries without visible residue
- Disinfection of electronic equipment in OTs, AKD units, ICU, NICU, Recovery Room, Medical Laboratories, etc.
- Excellent in critical areas like OTs, HIV endangered spaces and dental clinics
- Rapid and effective disinfection of skin flora
- Prevents intra-operative growth
- Unique de-fatting property
- Deep penetration, inhibits intra-operative bacterial growth in the incision folds of the skin
- Rapid on-set and good residual effect
- Excellent skin tolerance

Packaging:

- Provided with convenient hand sprayer with 2 nozzle options - jet spray and mist spray.



Areas of Application:

Equipment Disinfection:

- Electronic devices, imaging machines, respirators, ventilators, cardiac monitors, emergency trolleys, CT scanners, haemodialysis systems, operation tables, lights, etc.
- Emergency disinfection of operation theatres in between surgeries, ICU, ICCU, NICU, etc.

Skin preparation:

- For prepping and post closure of skin. Has excellent cleansing and de-fatting effects.

Directions for use:

Equipment Disinfection:

Set the spray nozzle at "Spray" position and keeping the nozzle around 25 cm away from the equipment surface, spray **Fumispray®-Rapid** liberally. Allow the liquid to stay on the surface for 30 seconds and wipe away with a clean, lint-free cloth.

Skin Preparation:

About 10 minutes prior to incision, spray the operative and peri-operative area with **Fumispray®-Rapid** and paint for 2 minutes with sterile swabs (one swab for 2-3 bands) Allow the area to dry (about 30 seconds) Spray and re-paint the entire area as before and allow to dry Cover the disinfected area with sterile towels, leaving the incision site open



Download
Test Reports



Download
MSDS



PSK PHARMA PRIVATE LIMITED

121/22 Mahajenahalli, Shimoga Road, Harihar - 577601, Karnataka, India

ISO 13485 Medical Devices
Certified QMS

Tel: +91 96060 69014 | e-mail: weprotect@pskpharma.com | www.pskpharma.com



WHO-GMP Certified