


FUMISPRAY® - ANTI-FOG

Anti-fogging & Disinfecting Spray

**ANTI-FOG ACTION
PREVENTS MISTING & FOGGING**



Composition:

Each 100 g Contains:

Ethanol IP	10.0 g
2-Propanol IP	9.0 g
1- Propanol	6.0 g
Anti-fog agent & Surfactant	



Areas of Applications:

- To prevent fogging and misting on lenses & other clear transparent surfaces during use
- Endoscope lenses
- Spectacles
- Face shields

Directions for Use:

- Spray liberally onto sterile gauze pad and rub onto lens or transparent surface of equipment until the lens / surface is completely wetted with the liquid.
- Do not wipe off liquid & allow liquid to evaporate completely before using equipment

Precautions:

- Do not spray onto hot surfaces
- Keep in a cool, dark place away from sources of heat & light

Packaging:

- 100 ml Mist Spray Bottle



Home