

HEXODENT[®]

Chlorhexidine Mouthwash. Antiseptic Antiplaque



Composition:

Chlorhexidine Gluconate Solution IP diluted to Chlorhexidine Gluconate 0.2 % w/v in mint flavoured aqueous base.

Colours: Tartrazine & Brilliant Blue FCF

Applications:

- Gingivitis, Aphthous Ulceration & Candidal Infection.
- Periodontal Surgery
- Control of Dental Caries
- Maintenance of Oral Hygiene.

Salient Features:

- Highly effective in control and prevention of dental caries
- Controls gingival inflammation and bleeding
- Effective against wide range of micro-organisms including gram-positive and gram-negative bacteria
- Significantly reduces bacterial count in 5 minutes after rinsing
- Reduces extent of tooth decay

Packaging:

- 100 ML
- 100 ML with Box
- 1 Litre



Indications:

- Gingivitis
- Stomatitis
- Tonsillitis
- Pharyngitis
- Laryngitis
- Halitosis

Dilution:

- To be used at 50% dilution with water or undiluted based on duration of use
- 50% dilution: For prolonged use
- Undiluted: For normal use

Directions for Use:

- One cupful of Hexodent to be swished around in the mouth for one minute twice a day
- Do not rinse with water, eat or drink thereafter for 10 minutes
- For prolonged use, dilute Hexodent with an equal volume of water (50% dilution)
- Recommended only under professional guidance



CERT. NO.: 99-100-14857
ISO 9001:2008



PSK PHARMA PRIVATE LIMITED

121/22, Mahajenahalli, Shimoga Road, Harihar-577 601, Karnataka, India
Tel : +91 8192 24 0317, +91 8192 24 0318, Fax : +91 8192 24 0002

email : weprotect@pskpharma.com
website : www.pskpharma.com



GMP Certified